



W&H adaptor set for cleaning units and disinfectors

Tips for preparation in cleaning units and disinfectors

Cleaning adaptor for ultrasonic scaler tips

(prophylaxis, periodontology, endodontics, restorative dentistry & prosthetics)

- > All scaler tips with a M3 x 0.5 mm thread (internal thread – viewed from the tip) fit on the cleaning adaptor with the REF 06895000
- > EMS tips with a M3 x 0.5 mm thread (internal thread) e.g., EMS C
[EMS tips do not fit the quick-clamping system.]
- > NSK tips with a M3 x 0.5 mm thread (internal thread), compatible with EMS handpiece
All other NSK tips do not fit
- > Mectron tips with a M3 x 0.5 mm thread (internal thread), e.g., Mectron S7
- > Satelec tips are not compatible (M3 x 0.6 mm thread)

Cleaning adaptors for ultrasonic scaler tips (surgery)

- > Only W&H Piezomed instruments fit the adaptor set for cleaning units and disinfectors
- > Tips like Sirona SIROSON, NSK SG3, Satelec CS3 and Mectron OT 7 are not compatible due to the different thread sizes

Cleaning adaptor for sonic scaler tips

(prophylaxis, periodontology, endodontics, restorative dentistry & prosthetics, surgery)

- > All scaler tips with a M3 x 0.5 mm thread (external thread – viewed from the tip) fit on the cleaning adaptor with the REF 06894700
- > Sirona tips with a M3 x 0.5 mm thread (external thread) like SIROAIR
- > KaVo tips with a M3 x 0.5 mm thread (external thread) like KaVo tip 36/10
[KaVo tips do not fit the quick-clamping system.]
- > NSK tips with a M3 x 0.5 mm thread (external thread) like NSK S1
- > Brasseler/Komet tips with a M3 x 0.5 mm thread (external thread) like SF8878KD

External coolant tubes, spray clips, Y manifolds, etc., can also be cleaned

Errors and omissions excepted