



Correlation of irrigation tubings to W&H surgical units

implantmed

SI-915 / SI-923



REF 16929001 / 16929000




REF 04708900 / 04801300 / 04691400 / 04801200

implantMED


SI-95



135°C
SSS




04013900 (1 pcs)
Pump tubing with nipple




04363600 (6 pcs, 2.2 m)
Disposable irrigation tubing set

135°C
SSS




04014000 (3 pcs, 85 mm)
Pump tubing

2




00929300 (10 pcs, 2.2 m)
Disposable irrigation tubing

135°C
SSS




06338400
Irrigation spike with clamp

135°C
SSS




04013900 (1 pcs)
Pump tubing with nipple




04363600 (6 pcs, 2.2 m)
04364100 (6 pcs, 3.8 m)
Disposable irrigation tubing set

135°C
SSS




04014000 (3 pcs, 85 mm)
Pump tubing

2



00929300 (10 pcs, 2.2 m)
00929400 (10 pcs, 3.8 m)
Disposable irrigation tubing

135°C
SSS



06338400
Irrigation spike with clamp



Correlation of irrigation tubings to W&H surgical units

elcomed

SA-310




elcoMED

SA-200




135°C
SSS




04013900 (1 pcs)
Pump tubing with nipple

135°C
SSS




04014000 (3 pcs, 85 mm)
Pump tubing

②




00929300 (10 pcs, 2.2 m)
00929400 (10 pcs, 3.8 m)
Disposable irrigation tubing




04363600 (6 pcs, 2.2 m)
04364100 (6 pcs, 3.8 m)
Disposable irrigation tubing set

135°C
SSS




06338400
Irrigation spike with clamp

135°C
SSS




04035100 (1 pcs)
Pump tubing with nipple

135°C
SSS




04035200 (3 pcs, 190 mm)
Pump tubing

②




00929300 (10 pcs, 2.2 m)
00929400 (10 pcs, 3.8 m)
Disposable irrigation tubing

②




04366300 (4 pcs, 2.2 m)
04364800 (4 pcs, 3.8 m)
Disposable Y-Switch

②



04365300 (6 pcs, 2.2 m)
04365100 (6 pcs, 3.8 m)
Disposable irrigation tubing set

135°C
SSS



06338400
Irrigation spike with clamp