

# Guide des références SmartPeg

Édition Octobre 2018 | 25016-66

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg	
<b>A.B. DENTAL DEVICES</b>	Screw Type, Conical Groovy				3.75	100440	Type 32	
	Screw Type, Conical Groovy				4.5	100431	Type 27	
<b>ADDENTA</b>	Bifasico				4.0	100374	Type 14	
<b>ADIN DENTAL IMPLANT SYSTEM</b>	One				3.3	100390	Type 17	
	One				4.2	100390	Type 17	
	Swell				3.3	100480	Type 49	
	Swell				3.75	100480	Type 49	
	Swell				4.2	100480	Type 49	
	Swell				5.0	100480	Type 49	
	Swell				6.0	100480	Type 49	
	Touareg CloseFit				NP	100545	Type 62	
	Touareg CloseFit				RP	100543	Type 60	
	Touareg CloseFit				WP	100480	Type 49	
	Touareg-S				4.2	100480	Type 49	
	Touareg-S				5.0	100480	Type 49	
<b>AGS MEDIKAL</b>	Implance Dental Implant System, Bone Level					100378	Type 06	
	Implance Dental Implant System, Mini Type					100455	Type 38	
	Implance Dental Implant System, Tissue Level					100378	Type 06	
<b>ALLIANCE</b>	SBN		NP		3.0	100543	Type 60	New
	SBS		SP		5.0	100544	Type 61	New
<b>ALPHA BIO</b>	ATI					100425	Type 26	
	ATID				3.75	100425	Type 26	
	Neo CS				3.75	100583	Type 78	
	NICE 3.2 Conical connection				3.2	100554	Type 63	
	SPI					100425	Type 26	

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>ALPHA-DENT</b>	Active 3D				3.25	100480	Type 49
	Active 3D				3.75	100480	Type 49
	Active 3D				4.2	100480	Type 49
	Active 3D				5.0	100480	Type 49
	Active Conus				3.25	100480	Type 49
	Active Conus				3.75	100480	Type 49
	Active Conus				4.2	100480	Type 49
	Active Conus				5.0	100480	Type 49
	Active Plus				4.2	100480	Type 49
	Active Plus				5.0	100480	Type 49
	Classic 3D				3.3	100480	Type 49
	Classic 3D				3.75	100480	Type 49
	Classic 3D				4.2	100480	Type 49
	Classic 3D				5.0	100480	Type 49
	Classic Conus				3.3	100480	Type 49
	Classic Conus				3.75	100480	Type 49
	Classic Conus				4.2	100480	Type 49
	Classic Conus				5.0	100480	Type 49
<b>ALPHA-GATE</b>	Bioactive				3.75	100446	Type 35
<b>ALPHATECH</b>					3.4	100477	Type 46
					3.8	100477	Type 46
					4.3	100477	Type 46
					5.0	100477	Type 46
					6.0	100477	Type 46
					7.0	100477	Type 46
<b>ANKYLOS</b>	C/X		A			100388	Type 16
	C/X		B			100388	Type 16
	C/X		C			100388	Type 16
<b>ANTHOGYR</b>	Anthofit External Hexagon HE					100353	Type 01
	Anthofit Internal Octagon OI					100398	Type 18
	Axiom PX				3.4	100455	Type 38

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>ANTHOGYR</b>	Axiom PX			4.0		100455	Type 38
	Axiom PX			4.6		100455	Type 38
	Axiom PX			5.2		100455	Type 38
	Axiom REG			3.4		100455	Type 38
	Axiom REG			4.0		100455	Type 38
	Axiom REG			4.6		100455	Type 38
	Axiom REG			5.2		100455	Type 38
	Ossfit					100350	Type 04
<b>ARDS IMPLANTS</b>	Classic			All		100431	Type 27
	Premium			All		100431	Type 27
	Smart			All		100431	Type 27
<b>AVINENT IMPLANT SYSTEM</b>	Coral External			3.5		100436	Type 30
	Coral External			4.1		100353	Type 01
	Coral External			5.1		100357	Type 03
	Internal Hexagon Coral			3.5		100440	Type 32
	Internal Hexagon Coral			4.1		100440	Type 32
	Internal Hexagon Ocean			3.5		100440	Type 32
	Internal Hexagon Ocean			4.1		100440	Type 32
	Ocean External			3.5		100436	Type 30
	Ocean External			4.1		100353	Type 01
<b>AXIS</b>	Ceramic			4.0		100455	Type 38
<b>B&amp;B DENTAL</b>	Duravit SLIM Implant			3.0		100404	Type 22
	Duravit SLIM Implant			3.4		100404	Type 22
	Duravit 3P Implant			3.5		100425	Type 26
	Duravit 3P Implant			4.0		100425	Type 26
	Duravit 3P Implant			4.5		100425	Type 26
	Duravit 3P Implant			5.0		100425	Type 26
	Duravit EV Implant			4.5		100425	Type 26
	Duravit EV Implant			4.0		100425	Type 26
	Duravit EV Implant			4.5		100425	Type 26
	Duravit EV Implant			5.0		100425	Type 26

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>B&amp;B DENTAL</b>	Duravit WIDE Implant			5.0		100425	Type 26
	Duravit WIDE Implant			6.0		100425	Type 26
	DURAVIT WIDE Implant			6.5		100425	Type 26
	DURAVIT WIDE Implant			7.0		100425	Type 26
<b>BEGO IMPLANT SYSTEMS</b>	Semados Mini			2.7		100355	Type 02
	Semados Mini			2.9		100355	Type 02
	Semados Mini			3.1		100355	Type 02
	Semados RI-line			3.75		100425	Type 26
	Semados RI-line			4.1		100425	Type 26
	Semados RI-line			4.5		100425	Type 26
	Semados RI-line			5.5		100425	Type 26
	Semados RS/RSX			3.0		100481	Type 50
	Semados RS/RSX			3.75		100425	Type 26
	Semados RS/RSX			4.1		100425	Type 26
	Semados RS/RSX			4.5		100425	Type 26
	Semados RS/RSX			5.5		100425	Type 26
	Semados S-line			3.25		100425	Type 26
	Semados S-line			3.75		100425	Type 26
	Semados S-line			4.1		100425	Type 26
	Semados S-line			4.5		100425	Type 26
Semados S-line			5.5		100425	Type 26	
<b>BIO3</b>	Implant Advanced ATI			3.3		100480	Type 49
	Implant Advanced ATI			3.8		100480	Type 49
	Implant Advanced ATI			4.2		100480	Type 49
	Implant Advanced ATI			5.0		100480	Type 49
	Implant Progressive PTI			3.3		100480	Type 49
	Implant Progressive PTI			3.8		100480	Type 49
	Implant Progressive PTI			4.2		100480	Type 49
	Implant Progressive PTI			5.0		100480	Type 49
<b>BIODENTA</b>	Tissue level		NP			100402	Type 21
	Tissue level		RP			100478	Type 47
	Tissue level		WP			100478	Type 47

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>BIOHORIZONS</b>	External				3.5	100436	Type 30
	External				4.0	100353	Type 01
	External				5.0	100353	Type 01
	External				6.0	100353	Type 01
	Internal		Yellow			100440	Type 32
	Internal		Green			100431	Type 27
	Internal		Blue			100442	Type 33
	Internal Laser Lok				3.0	100425	Type 26
	Single Stage		Yellow			100440	Type 32
	Single Stage		Green			100431	Type 27
	Single Stage		Blue			100442	Type 33
	Tapered Internal		Yellow			100440	Type 32
	Tapered Internal		Green			100431	Type 27
	Tapered Internal		Blue			100442	Type 33
<b>BIOINFINITY</b>	Junior Plateforme				3.2	100455	Type 38
	Standard Plateforme				4.8	100425	Type 26
<b>BIOMET 3i</b>	OSSEOTITE Certain		XP		4/5	100386	Type 15
	OSSEOTITE Certain		XP		5/6	100386	Type 15
	OSSEOTITE Tapered Certain		Purple			100472	Type 45
	OSSEOTITE Tapered Certain		Blue			100386	Type 15
	NanoTite Certain Microminiplant				3.25	100472	Type 45
	NanoTite External Hex		Purple			100353	Type 01
	NanoTite External Hex		Blue			100353	Type 01
	NanoTite External Hex		Yellow			100353	Type 01
	NanoTite External Hex		Green			100353	Type 01
	NanoTite External Hex				XP 3/4	100353	Type 01
	NanoTite External Hex				XP 4/5	100353	Type 01
	NanoTite External Hex				XP 5/6	100353	Type 01
	NanoTite External Hex MicroMiniplant		Purple			100353	Type 01
	NanoTite Tapered Certain		Purple			100472	Type 45
	NanoTite Tapered Certain		Blue			100386	Type 15
	NanoTite Tapered Certain		Yellow			100386	Type 15
	NanoTite Tapered Certain		Green			100386	Type 15

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>BIOMET 3i</b>	NanoTite Tapered Certain Prevail		Purple			100472	Type 45
	NanoTite Tapered Certain Prevail		Blue			100386	Type 15
	NanoTite Tapered Certain Prevail		Yellow			100386	Type 15
	OSSEOTITE Certain Microminiplant			3.25		100472	Type 45
	OSSEOTITE External Hex		Purple			100353	Type 01
	OSSEOTITE External Hex		Blue			100353	Type 01
	OSSEOTITE External Hex		Yellow			100353	Type 01
	OSSEOTITE External Hex		Green			100353	Type 01
	OSSEOTITE External Hex			XP 3/4		100353	Type 01
	OSSEOTITE External Hex			XP 4/5		100353	Type 01
	OSSEOTITE External Hex			XP 5/6		100353	Type 01
	OSSEOTITE External Hex MicroMiniplant		Purple			100353	Type 01
	OSSEOTITE Tapered Certain		Yellow			100386	Type 15
	OSSEOTITE Tapered Certain		Green			100386	Type 15
	OSSEOTITE Tapered Certain Prevail		Purple			100472	Type 45
	OSSEOTITE Tapered Certain Prevail		Blue			100386	Type 15
	OSSEOTITE Tapered Certain Prevail		Yellow			100386	Type 15
	T3 Non-Plateforme Switched		Purple			100472	Type 45
	T3 Non-Plateforme Switched		Blue			100386	Type 15
	T3 Non-Plateforme Switched		Yellow			100386	Type 15
	T3 Non-Plateforme Switched		Green			100386	Type 15
	T3 Plateforme Switched		Purple			100472	Type 45
	T3 Plateforme Switched		Blue			100386	Type 15
	T3 Plateforme Switched		Yellow			100386	Type 15
<b>BIOTEC SRL</b>	EN				All	100436	Type 30
	ER				All	100353	Type 01
	IM				All	100431	Type 27
	IR				All	100440	Type 32
	KR				3.0	100455	Type 38
	KW				3.6	100440	Type 32
	SR				All	100378	Type 06
	TWO				All	100440	Type 32
<b>BIOTECH DENTAL</b>	Contact K30				3.0	100581	Type 75

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>BIOTECH DENTAL</b>	Kontakt				3.6	100582	Type 76
	Kontakt				4.2	100582	Type 76
	Kontakt				4.8	100582	Type 76
	Kontakt				5.4	100582	Type 76
	Weego.FR				3.6	100583	Type 78
	Weego.FR				4.2	100583	Type 78
	Weego.FR				4.8	100583	Type 78
<b>BONE SYSTEM</b>	2P				5.0	100447	Type 36
<b>BRAINBASE</b>			C			100478	Type 47
	Mytis Arrow Implants		B			100350	Type 04
	Mytis Arrow Implants		BW			100350	Type 04
	Mytis Arrow Implants		E			100350	Type 04
	Mytis Arrow Implants		EW			100350	Type 04
<b>BREDDENT</b>	blueSKY					100480	Type 49
<b>BTI BIOTECHNOLOGY INSTITUTE</b>	Externa		Universal			100353	Type 01
	Externa		Universal plus			100353	Type 01
	Interna		Universal			100431	Type 27
	Interna		Universal plus			100431	Type 27
	Interna		Wide			100431	Type 27
	Tiny			2.5		100368	Type 09
	Tiny			3.75		100368	Type 09
<b>BTLOCK</b>	CV1				3.75	100480	Type 49
	CV1				4.5	100480	Type 49
	CV1				5.5	100480	Type 49
	CV1				6.5	100480	Type 49
	CV2				3.75	100480	Type 49
	CV2				4.5	100480	Type 49
	CV2				5.5	100480	Type 49
	CV2				6.5	100480	Type 49
	CV3				3.75	100480	Type 49

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>BTLOCK</b>	CV3			4.5		100480	Type 49
	CV3			5.5		100480	Type 49
	CV3			6.5		100480	Type 49
<b>CAMLOG</b>	Camlog Cylinder-line			3.3		100432	Type 28
	Camlog Cylinder-line			3.8		100411	Type 23
	Camlog Cylinder-line			4.3		100411	Type 23
	Camlog Cylinder-line			5.0		100412	Type 24
	Camlog Cylinder-line			6.0		100412	Type 24
	Camlog Root-line Promote			3.8		100411	Type 23
	Camlog Root-line Promote			4.3		100411	Type 23
	Camlog Root-line Promote			5.0		100412	Type 24
	Camlog Root-line Promote			6.0		100412	Type 24
	Camlog Screw-line			3.3		100432	Type 28
	Camlog Screw-line			3.8		100411	Type 23
	Camlog Screw-line			4.3		100411	Type 23
	Camlog Screw-line			5.0		100412	Type 24
	Camlog Screw-line			6.0		100412	Type 24
	Conelog Screw-line			3.3		100557	Type 65
	Conelog Screw-line			3.8		100541	Type 58
	Conelog Screw-line			4.3		100541	Type 58
	Conelog Screw-line			5.0		100380	Type 07
	iSy (direct to implant)			All		100497	Type 57
	iSy (to implant base)			All		100446	Type 35
<b>CHANG GUNG</b>	ComMed			3.5		100481	Type 50
	ComMed			4.0		100481	Type 50
	ComMed			4.5		100380	Type 07
	ComMed			5.0		100380	Type 07
	ComMed			5.5		100380	Type 07
	ComMed			6.0		100380	Type 07
<b>CHAORUM</b>	External			3.3		100353	Type 01
	Internal			4.0		100478	Type 47
	Internal			4.8		100478	Type 47



Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>COCHLEAR</b>	Baha & Vistafix Flange Fixture Hexagonal Connection				3.0	100368	Type 09
	Baha & Vistafix Flange Fixture Hexagonal Connection				4.0	100368	Type 09
	Baha BI 300				3.0	100436	Type 30
	Baha BI 300				4.0	100436	Type 30
	Baha BIA 300				6.0	100487	Type 55
	Baha BIA 300				9.0	100487	Type 55
	Baha BIA 400				6.0	100487	Type 55
	Baha BIA 400				8.0	100487	Type 55
	Baha BIA 400				10.0	100487	Type 55
	Baha BIA 400				12.0	100487	Type 55
	Baha BIA 400				14.0	100487	Type 55
	Vistafix VXI 300				3.0	100436	Type 30
	Vistafix VXI 300				4.0	100436	Type 30
<b>COMPRESS</b>					4.0	100355	Type 02
<b>CORTEX</b>	Classix				3.3	100425	Type 26
	Classix				3.8	100425	Type 26
	Classix				4.2	100425	Type 26
	Classix				5.0	100425	Type 26
	Classix				6.0	100425	Type 26
	Dynamix				3.3	100425	Type 26
	Dynamix				3.8	100425	Type 26
	Dynamix				4.2	100425	Type 26
	Dynamix				5.0	100425	Type 26
	Dynamix				6.0	100425	Type 26
	Saturn				3.3	100425	Type 26
	Saturn				3.8	100425	Type 26
	Saturn				4.2	100425	Type 26
	Saturn				5.0	100425	Type 26
Saturn				6.0	100425	Type 26	
<b>COWELLMEDI</b>	Atlas External				3.3	100355	Type 02
	Atlas External				4.1	100353	Type 01

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>C-TECH</b>	Atlas External				5.0	100353	Type 01
	Atlas Submerged				4.0	100402	Type 21
	EL				3.5	100425	Type 26
	EL				4.3	100425	Type 26
	EL				5.1	100425	Type 26
<b>DENTAL TECH</b>			WP			100478	Type 47
<b>DENTAURUM</b>	Tiologic		S			100472	Type 45
	Tiologic		M			100472	Type 45
	Tiologic		L			100472	Type 45
<b>DENTEGRIS</b>	Straight & Tapered				3.3	100425	Type 26
	Straight & Tapered				3.75	100440	Type 32
	Straight & Tapered				4.5	100440	Type 32
	Straight & Tapered				4.75	100440	Type 32
	Straight & Tapered				5.5	100478	Type 47
<b>DENTIS DENTAL</b>	Internal				3.7	100350	Type 04
	Internal				4.3	100350	Type 04
	Internal				4.8	100350	Type 04
	OneQ S-clean Narrow (SLA surface/Submerged)				3.0	100404	Type 22
	OneQ S-clean Narrow (SLA surface/Submerged)				3.3	100404	Type 22
	OneQ S-clean Regular (SLA surface/Submerged)				3.7	100372	Type 11
	OneQ S-clean Regular (SLA surface/Submerged)				4.2	100372	Type 11
	OneQ S-clean Regular (SLA surface/Submerged)				4.7	100372	Type 11
	OneQ S-clean Regular (SLA surface/Submerged)				5.2	100372	Type 11
	OneQ S-clean Wide (SLA surface/Submerged)				6.0	100372	Type 11
	OneQ S-clean Wide (SLA surface/Submerged)				7.0	100372	Type 11
	OneQ S-clean Wide (SLA surface/Submerged)				8.0	100372	Type 11
	Save Submerged				5.5	100380	Type 07
	Save Submerged				6.0	100380	Type 07
	S-clean Mini (Submerged)				3.7	100372	Type 11
	S-clean Regular (Submerged)				4.1	100372	Type 11
	S-clean Regular (Submerged)				4.3	100372	Type 11

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>DENTIS DENTAL</b>	S-clean Wide (Submerged)			4.8		100372	Type 11
<b>DENTIUM</b>	Genoss			3.0		100554	Type 63
	Genoss			3.6		100402	Type 21
	Genoss			4.3		100402	Type 21
	Genoss			5.0		100402	Type 21
	Implantium		Yellow 3.6	3.4		100378	Type 06
	Implantium		Green 4.0	3.8		100378	Type 06
	Implantium		Blue 4.5	4.3		100378	Type 06
	Simple Line			3.8		100380	Type 07
	Simple Line			4.3		100380	Type 07
	Simple Line			4.8		100380	Type 07
	Super Line			4.0		100380	Type 07
	Super Line			7.0		100380	Type 07
<b>DENTSPLY IMPLANTS</b>	ANKYLOS C/X Implants		A			100388	Type 16
	ANKYLOS C/X Implants		B			100388	Type 16
	ANKYLOS C/X Implants		C			100388	Type 16
	ASTRA TECH Implant system MicroThread			3.5		100376	Type 05
	ASTRA TECH Implant system MicroThread			4.0		100378	Type 06
	ASTRA TECH Implant system MicroThread			4.5		100380	Type 07
	ASTRA TECH Implant system MicroThread			5.0		100380	Type 07
	ASTRA TECH Implant system OsseoSpeed Aqua			3.5 S		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed Aqua			4.0 S		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25212	Green	3.0 S	8	100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25213	Green	3.0 S	9	100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25214	Green	3.0 S	11	100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25215	Green	3.0 S	13	100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25216	Green	3.0 S	15	100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25221	Purple	3.6 S	6	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25222	Purple	3.6 S	8	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25223	Purple	3.6 S	9	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25224	Purple	3.6 S	11	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25225	Purple	3.6 S	13	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25226	Purple	3.6 S	15	100455	Type 38

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>DENTSPLY IMPLANTS</b>	ASTRA TECH Implant system OsseoSpeed EV	25227	Purple	3.6 S	17	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25231	Yellow	4.2 S	6	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25232	Yellow	4.2 S	8	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25233	Yellow	4.2 S	9	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25234	Yellow	4.2 S	11	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25235	Yellow	4.2 S	13	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25236	Yellow	4.2 S	15	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25237	Yellow	4.2 S	17	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25241	Blue	4.8 S	6	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25242	Blue	4.8 S	8	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25243	Blue	4.8 S	9	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25244	Blue	4.8 S	11	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25245	Blue	4.8 S	13	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25246	Blue	4.8 S	15	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25247	Blue	4.8 S	17	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25251	Brown	5.4 S	6	100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25252	Brown	5.4 S	8	100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25253	Brown	5.4 S	9	100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25254	Brown	5.4 S	11	100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25255	Brown	5.4 S	13	100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25256	Brown	5.4 S	15	100462	Type 40
	ASTRA TECH Implant system OsseoSpeed EV	25262	Yellow	4.2 C	8	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25263	Yellow	4.2 C	9	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25264	Yellow	4.2 C	11	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25265	Yellow	4.2 C	13	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25266	Yellow	4.2 C	15	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25267	Yellow	4.2 C	17	100480	Type 49
	ASTRA TECH Implant system OsseoSpeed EV	25272	Blue	4.8 C	8	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25273	Blue	4.8 C	9	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25274	Blue	4.8 C	11	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25275	Blue	4.8 C	13	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25276	Blue	4.8 C	15	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed EV	25277	Blue	4.8 C	17	100478	Type 47
	ASTRA TECH Implant system OsseoSpeed Lilac			4.5		100380	Type 07
	ASTRA TECH Implant system OsseoSpeed Lilac			5.0		100380	Type 07

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>DENTSPLY IMPLANTS</b>	ASTRA TECH Implant system OsseoSpeed Lilac			5.0S		100380	Type 07
	ASTRA TECH Implant system OsseoSpeed Profile EV	25452	Yellow	4.2 PC	8	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25453	Yellow	4.2 PC	9	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25454	Yellow	4.2 PC	11	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25455	Yellow	4.2 PC	13	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25456	Yellow	4.2 PC	15	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25457	Yellow	4.2 PC	17	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25458	Blue	4.8 PC	8	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25459	Blue	4.8 PC	9	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25460	Blue	4.8 PC	11	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25461	Blue	4.8 PC	13	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25462	Blue	4.8 PC	15	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25463	Blue	4.8 PC	17	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25440	Yellow	4.2 PS	8	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25441	Yellow	4.2 PS	9	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25442	Yellow	4.2 PS	11	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25443	Yellow	4.2 PS	13	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25444	Yellow	4.2 PS	15	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25445	Yellow	4.2 PS	17	100370	Type 10
	ASTRA TECH Implant system OsseoSpeed Profile EV	25446	Blue	4.8 PS	8	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25447	Blue	4.8 PS	9	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25448	Blue	4.8 PS	11	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25449	Blue	4.8 PS	13	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25450	Blue	4.8 PS	15	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed Profile EV	25451	Blue	4.8 PS	17	100440	Type 32
	ASTRA TECH Implant system OsseoSpeed TX Aqua			3.5S		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed TX Aqua			4.0S		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed TX Lilac			4.5		100380	Type 07
	ASTRA TECH Implant system OsseoSpeed TX Lilac			5.0		100380	Type 07
	ASTRA TECH Implant system OsseoSpeed TX Lilac			5.0 S		100380	Type 07
	ASTRA TECH Implant system OsseoSpeed TX Profile			4.5		100378	Type 06
	ASTRA TECH Implant system OsseoSpeed TX Profile			5.0		100378	Type 06
	ASTRA TECH Implant system OsseoSpeed TX Profile			5.0 S		100378	Type 06
	ASTRA TECH Implant system OsseoSpeed TX Yellow			3.0 S		100481	Type 50
	ASTRA TECH Implant system OsseoSpeed Yellow			3.0 S		100481	Type 50

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg	
<b>DENTSPLY IMPLANTS</b>	XiVe		TG			100366	Type 08	
	XiVe				3.0	100404	Type 22	
	XiVe				3.4	100402	Type 21	
	XiVe				3.8	100472	Type 45	
	XiVe				4.5	100472	Type 45	
	XiVe				5.5	100472	Type 45	
<b>DIO</b>	Bone Level UF (II) Narrow		1.7 Hex			100404	Type 22	New
	Bone Level UF (II) Regular		2.5 Hex			100378	Type 06	New
	Bone Level UF (II) Wide		2.5 Hex			100378	Type 06	New
	External FSN				3.3	100355	Type 02	
	External FSN				3.75	100353	Type 01	
	External FTN Narrow		2.4 Hex			100355	Type 02	*
	External FTN Regular		2.7 Hex			100353	Type 01	*
	External FTN T-Wide		2.7 Hex			100353	Type 01	*
	External FTN Wide		3.4 Hex			100357	Type 03	*
	SM Submerged				3.8	100480	Type 49	
	SM Submerged				4.5	100478	Type 47	
	SM Submerged				5.3	100478	Type 47	
	Tissue Level UF (II) P4.8		2.9 Octa			100350	Type 04	New
	Tissue Level UF (II) P6.0		2.9 Octa			100350	Type 04	New
<b>DTI</b>					3.5	100497	Type 57	
					4.2	100380	Type 07	
<b>DYNA DENTAL</b>	Dyna Helix DC				3.2	100431	Type 27	
	Dyna Helix DC				3.6	100431	Type 27	
	Dyna Helix DC				4.2	100431	Type 27	
	Dyna Helix DC				5.0	100431	Type 27	
	Dyna Helix ST				3.6	100431	Type 27	
	Dyna Helix ST				4.2	100431	Type 27	
	Dyna Helix ST				5.0	100431	Type 27	
<b>ECKERMANN</b>	All Spiral					100390	Type 17	
	Hexagono Ext Hex					100353	Type 01	

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>ECKERMANN</b>	Tripto					100480	Type 49
<b>EMFILS</b>	Colosso IH ECI Evolution			3.3		100372	Type 11
	Colosso IH ECI Evolution			4.0		100372	Type 11
	Colosso IH ECI Evolution			5.0		100372	Type 11
	ICM-4010			4.0		100484	Type 53
	PAT5-4010			4.0		100484	Type 53
<b>ENTIFIC</b>	Flange Fixture					100353	Type 01
<b>EUROTEKNIKA</b>	Aesthetica+2		3.6	4.2		100544	Type 61
	Aesthetica+2		3.6	4.8		100544	Type 61
	Aesthetica+2		4.1	4.2		100544	Type 61
	Aesthetica+2		4.1	4.8		100544	Type 61
	Aesthetica+2		4.8	4.8		100544	Type 61
	Aesthetica+2		4.8	6.5		100544	Type 61
	Naturall			3.5		100455	Type 38
	Naturall			4.0		100378	Type 06
	Naturall			4.5		100378	Type 06
	Naturall			5.0		100378	Type 06
	Uneva+			3.6		100353	Type 01
	Uneva+			4.1		100353	Type 01
	Uneva+			4.8		100353	Type 01
	Uneva+			6.0		100353	Type 01
<b>FAIR IMPLANT</b>	FairTwo			3.5		100484	Type 53
	FairTwo			4.2		100484	Type 53
	FairTwo			5.0		100378	Type 06
	FairTwo			6.0		100378	Type 06
<b>GC IMPLANT</b>	Aadva		Narrow			100455	Type 38
	Aadva		Regular			100541	Type 58
	Aadva		Wide			100541	Type 58
	External Plus			3.0		100402	Type 21
	External Plus			3.4		100402	Type 21

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>GC IMPLANT</b>	External Plus				3.8	100353	Type 01
	External Plus				5.0	100353	Type 01
	GENESIO					100402	Type 21
	Internal Plus				3.0	100402	Type 21
	Internal Plus				3.4	100402	Type 21
	SETIO					100353	Type 01
<b>GHIMAS</b>	PFON EVO SERIES				3.6	100455	Type 38
	PFON EVO SERIES				4.0	100455	Type 38
	PFON EVO SERIES				5.0	100380	Type 07
<b>HIOSSEN</b>	ETIII		Mini		3.5	100376	Type 05
	ETIII				4.0	100478	Type 47
	ETIII				4.5	100478	Type 47
	ETIII				5.0	100478	Type 47
	ETIII		UW		6.0	100378	Type 06
	ETIII		UW		7.0	100378	Type 06
	HGII		Mini			100376	Type 05
	HGII		Regular			100378	Type 06
	HGII		Wide			100378	Type 06
	HGII		Ultra Wide			100378	Type 06
	HGIII		Mini			100376	Type 05
	HSII		Regular			100350	Type 04
	HSII		Wide			100350	Type 04
<b>IBS INNO BIO SURG</b>	Magicore				5.5	100455	Type 38
<b>IDI EVOLUTION</b>	Vite				3.75	100398	Type 18
	Vite				4.0	100398	Type 18
<b>IMMEDIATELOAD</b>	IM40 Immediate Conico				4.0	100480	Type 49
<b>IMPLANT DIRECT</b>	Legacy1		Blue			100455	Type 38
	Legacy1		Green			100440	Type 32
	Legacy1		Purple			100431	Type 27



Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>IMPLANT DIRECT</b>	Legacy1		Yellow			100442	Type 33
	Legacy2		Blue			100455	Type 38
	Legacy2		Green			100440	Type 32
	Legacy2		Purple			100431	Type 27
	Legacy2		Yellow			100442	Type 33
	Legacy3		Blue			100455	Type 38
	Legacy3		Green			100440	Type 32
	Legacy3		Purple			100431	Type 27
	Legacy3		Yellow			100442	Type 33
	Legacy4		Blue			100455	Type 38
	Legacy4		Green			100440	Type 32
	Legacy4		Purple			100431	Type 27
	Legacy4		Yellow			100442	Type 33
	ScrewPlant				3.7	100425	Type 26
	ScrewPlant				4.7	100431	Type 27
	ScrewPlant				5.7	100431	Type 27
<b>IMTEC</b>	Hexed-head				3.75	100353	Type 01
	Hexed-head				4.75	100353	Type 01
<b>INTRA-LOCK SYSTEM</b>	Internal		SQ			100378	Type 06
	Internal		Standard			100478	Type 47
	Internal		Wide			100478	Type 47
<b>IRES</b>	IMAX NHSE				3.3	100355	Type 02
	IMAX NHSE				3.7	100353	Type 01
	IMAX NHSE				4.1	100353	Type 01
	IMAX NHSE				5.2	100353	Type 01
	IMAX NHSI				3.7	100431	Type 27
	IMAX NHSI				4.1	100431	Type 27
	IMAX NHSI				5.2	100431	Type 27
	Shapeone S1B				3.1	100431	Type 27
	Shapeone S1B				4.1	100431	Type 27
	Shapeone S1B				4.7	100431	Type 27
	Shapeone S1EH				3.1	100353	Type 01

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>IRES</b>	Shapeone S1EH			4.1		100353	Type 01
	Shapeone S1EH			4.7		100353	Type 01
	Shapeone S1T			3.7		100350	Type 04
	Shapeone S1T			4.1		100350	Type 04
	Shapeone S1T			4.7		100350	Type 04
	Shapeone S1TN			3.7		100431	Type 27
	Shapeone S1TN			4.1		100431	Type 27
	Shapeone S1TN			4.7		100431	Type 27
	Volution SVB			3.7		100431	Type 27
	Volution SVB			4.1		100431	Type 27
	Volution SVB			4.7		100431	Type 27
	<b>ITERUM</b>	Conical			3.75		100425
<b>JDENTALCARE</b>	JDEvolution			3.2		100455	Type 38
	JDEvolution			3.7		100440	Type 32
	JDEvolution			4.3		100440	Type 32
	JDEvolution			5.0		100440	Type 32
<b>JMM</b>	POI 2 piece			3.2		100398	Type 18
	POI 2 piece			3.7		100398	Type 18
	POI 2 piece			4.2		100398	Type 18
	POIEX			3.7		100386	Type 15
	POIEX			4.2		100480	Type 49
	POIEX			4.7		100478	Type 47
	POIEX			5.2		100478	Type 47
<b>KENTEC</b>	SB1			3.3		100378	Type 06
	SB2			4.2		100378	Type 06
<b>KEYSTONE DENTAL</b>	Genesis			3.8		100431	Type 27
	Genesis			4.5		100442	Type 33
	Genesis			5.5		100442	Type 33
	Genesis			6.5		100442	Type 33
	Prima Connex		SD	3.5		100431	Type 27

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>KEYSTONE DENTAL</b>	Prima Connex		RD	4.1		100442	Type 33
	Prima Connex		WD	5.0		100442	Type 33
	Renova RDI			4.5		100431	Type 27
	Renova SDI			3.75		100425	Type 26
	Restore			3.3		100436	Type 30
	Restore			3.75		100353	Type 01
	Restore			4.0		100353	Type 01
	Restore			5.0		100357	Type 03
	Restore			6.0		100357	Type 03
<b>KLOCKNER IMPLANT SYSTEM</b>	MV VEGA Implant			3.0		100497	Type 57
	NV VEGA Implant			3.5		100497	Type 57
	RV VEGA Implant			4.0		100425	Type 26
	RV VEGA Implant			4.5		100425	Type 26
	Essential Cone			3.5		100480	Type 49
	Essential Cone			4.0		100480	Type 49
	Essential Cone			4.5		100480	Type 49
	Essential Cone			4.8		100480	Type 49
	KL Implant NP			3.3		100355	Type 02
	KL Implant NP			3.7		100353	Type 01
	KL Implant NP			4.2		100353	Type 01
	KL Implant WP			4.7		100357	Type 03
	S3M					100436	Type 30
	SK2					100436	Type 30
<b>KYOCERA</b>	FINESIA BL		NP			100545	Type 62
	FINESIA BL		RP			100431	Type 27
	FINESIA BL		WP			100350	Type 04
	FINESIA TL		RP			100350	Type 04
	FINESIA TL		WP			100350	Type 04
<b>LASAK</b>	BioniQ, QN (Q-Lock Narrow)					100404	Type 22
	BioniQ, QR (Q-Lock Regular)					100497	Type 57
	Impladent			2.9		100355	Type 02
	Impladent			3.7		100447	Type 36

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>LASAK</b>	Impladent			5.0		100448	Type 37
<b>LEADER</b>	Implus EH			3.3		100353	Type 01
	Implus EH			3.75		100353	Type 01
	Implus EH			5.0		100353	Type 01
	Implus IH			3.3		100425	Type 26
	Tixos IH			3.75		100425	Type 26
	Tixos IH			4.5		100425	Type 26
	Tixos IH			5.5		100425	Type 26
	Tixos ITMC			3.75		100497	Type 57
	Tixos ITMC			4.5		100497	Type 57
	Tixos ITMC			5.5		100497	Type 57
<b>MEDENTIKA</b>	M-Implant®			3.5		100455	Type 38
	M-Implant®			4.0		100455	Type 38
	M-Implant®			4.5		100455	Type 38
	M-Implant®			5.0		100455	Type 38
<b>MEDICAL INSTINCT</b>	Bone Trust			3.4		100466	Type 44
	Bone Trust			4.0		100466	Type 44
	Bone Trust			5.0		100466	Type 44
	Bone Trust Plus Cone			3.4		100466	Type 44
	Bone Trust Plus Cone			4.0		100466	Type 44
<b>MEDIGMA TECHNOLOGIES</b>	Fix 1					100390	Type 17
	Fix-A-Dent					100480	Type 49
	FixTite					100480	Type 49
	FixTite-S					100480	Type 49
<b>MEGAGEN</b>	AnyOne			3.5		100380	Type 07
	AnyOne			4.0		100380	Type 07
	AnyOne External			4.5		100350	Type 04
	AnyRidge					100431	Type 27
	EZ Plus Internal		Small			100455	Type 38
	EZ Plus Internal		Regular			100380	Type 07

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>MEGAGEN</b>	EZ Plus Internal		Wide			100380	Type 07
	EZ Plus Internal			4.0		100380	Type 07
	EZ Plus Internal			4.5		100380	Type 07
	EZ Plus Internal			5.0		100380	Type 07
	Mini Internal			3.0		100481	Type 50
	Resque External			7.0		100357	Type 03
<b>MICRODENT</b>	Duraplant 2.2			4.1		100370	Type 10
	Microdent system					100353	Type 01
	MK			5.6		100357	Type 03
	MS			3.5		100436	Type 30
	MST			3.5		100436	Type 30
	Trylogic			3.8		100431	Type 27
	Trylogic			4.3		100431	Type 27
	Universal			4.1		100353	Type 01
	Universal			5.1		100353	Type 01
<b>MIS</b>	C1 Conical Connection		NP			100455	Type 38
	C1 Conical Connection		SP			100480	Type 49
	C1 Conical Connection		WP			100584	Type 79
	V3 Conical Connection		NP			100455	Type 38
	V3 Conical Connection		SP			100480	Type 49
	Internal Hexagon (seven, biocom)			3.75		100440	Type 32
	Internal Octagon (mistral)			4.75		100350	Type 04
	Seven		Narrow	3.3		100554	Type 63
<b>MOZO-GRAU</b>	MG-Inhex		Standard			100455	Type 38
	MG-Inhex		Maxi			100380	Type 07
	MG-Osseous		Mini			100353	Type 01
	MG-Osseous		Standard			100353	Type 01
	MG-Osseous		Maxi			100357	Type 03
<b>ND IMPLANTS</b>	Niko			3.6		100431	Type 27
	Niko			4.6		100380	Type 07

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg	
<b>NEMRIS</b>	Aesthura Classic			3.75		100446	Type 35	
	Aesthura Classic			4.75		100446	Type 35	
<b>NEOBIOTECH</b>	IS Fixture (RBM Surface)			4.0		100376	Type 05	
	IS Fixture (RBM Surface)			4.5		100376	Type 05	
	IS Fixture (RBM Surface)			4.8		100376	Type 05	
	IS S-Wide Fixture (RBM Surface)			5.5		100376	Type 05	
	IS S-Wide Fixture (RBM Surface)			6.0		100376	Type 05	
	IS S-Wide Fixture (RBM Surface)			7.0		100376	Type 05	
	IS-II active Fixture (S.L.A. Surface)			3.5		100376	Type 05	
	IS-II active Fixture (S.L.A. Surface)			4.0		100376	Type 05	
	IS-II active Fixture (S.L.A. Surface)			4.5		100376	Type 05	
	IS-II active Fixture (S.L.A. Surface)			5.0		100376	Type 05	
	IS-II active S-Wide Fixture (S.L.A. Surface)			5.5		100376	Type 05	
	IS-II active S-Wide Fixture (S.L.A. Surface)			6.0		100376	Type 05	
	IS-II active S-Wide Fixture (S.L.A. Surface)			7.0		100376	Type 05	
	IS-II Fixture (RBM Surface)			3.5		100376	Type 05	
	IS-II Fixture (RBM Surface)			4.0		100376	Type 05	
	IS-II Fixture (RBM Surface)			4.5		100376	Type 05	
	IS-II Fixture (RBM Surface)			5.0		100376	Type 05	
	IS-II S-Wide Fixture (RBM Surface)			5.5		100376	Type 05	
	IS-II S-Wide Fixture (RBM Surface)			6.0		100376	Type 05	
	IS-II S-Wide Fixture (RBM Surface)			7.0		100376	Type 05	
	IS-III active Fixture			3.5		100376	Type 05	
	IS-III active Fixture			4.0		100376	Type 05	
	IS-III active Fixture			4.5		100376	Type 05	
	IS-III active Fixture			5.0		100376	Type 05	
	IS-III active S-Wide Fixture			5.5		100376	Type 05	
	IS-III active S-Wide Fixture			6.0		100376	Type 05	
	IS-III active S-Wide Fixture			7.0		100376	Type 05	
	Neoplant External			4.0		100353	Type 01	
	Neoplant External			5.0		100353	Type 01	
	<b>NEODENT</b>	Alvim CM	109.657		3.5	8	100388	Type 16
Alvim CM		109.658		3.5	10	100388	Type 16	*

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>NEODENT</b>	Alvim CM	109.659			3.5 11.5	100388	Type 16 *
	Alvim CM	109.621			3.5 13	100388	Type 16 *
	Alvim CM	109.660			3.5 16	100388	Type 16 *
	Alvim CM	109.647			4.3 8	100388	Type 16 *
	Alvim CM	109.648			4.3 10	100388	Type 16 *
	Alvim CM	109.649			4.3 11.5	100388	Type 16 *
	Alvim CM	109.622			4.3 13	100388	Type 16 *
	Alvim CM	109.651			4.3 16	100388	Type 16 *
	Alvim CM	109.652			5.0 8	100388	Type 16 *
	Alvim CM	109.653			5.0 10	100388	Type 16 *
	Alvim CM	109.654			5.0 11.5	100388	Type 16 *
	Alvim CM	109.655			5.0 13	100388	Type 16 *
	Alvim CM	109.656			5.0 16	100388	Type 16 *
	Alvim Ti	109.321	4.3		4.3 10	100353	Type 01 *
	Alvim Ti	109.322	4.3		4.3 13	100353	Type 01 *
	Alvim Ti	109.323	4.3		4.3 16	100353	Type 01 *
	Alvim Ti	109.324	5.0		5.0 10	100357	Type 03 *
	Alvim Ti	109.325	5.0		5.0 13	100357	Type 03 *
	Alvim Ti	109.326	5.0		5.0 16	100357	Type 03 *
	Alvim II Plus				3.75	100402	Type 21
	Alvim II Plus				4.3	100402	Type 21
	Alvim II Plus				5.0	100402	Type 21
	Drive CM Acqua	109.692			3.5 8	100479	Type 48 *
	Drive CM Acqua	109.682			3.5 10	100388	Type 16 *
	Drive CM Acqua	109.693			3.5 11.5	100479	Type 48 *
	Drive CM Acqua	109.683			3.5 13	100388	Type 16 *
	Drive CM Acqua	109.684			3.5 16	100388	Type 16 *
	Drive CM Acqua	109.896			3.5 18	100479	Type 48 *
	Drive CM Acqua	109.689			4.3 8	100479	Type 48 *
	Drive CM Acqua	109.688			4.3 10	100479	Type 48 *
	Drive CM Acqua	109.627			4.3 11.5	100479	Type 48 *
	Drive CM Acqua	109.628			4.3 13	100479	Type 48 *
	Drive CM Acqua	109.629			4.3 16	100479	Type 48 *
	Drive CM Acqua	109.897			4.3 18	100479	Type 48 *
	Drive CM Acqua	109.690			5.0 8	100479	Type 48 *

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg	
<b>NEODENT</b>	Drive CM Acqua	109.685		5.0	10	100479	Type 48	*
	Drive CM Acqua	109.691		5.0	11.5	100479	Type 48	*
	Drive CM Acqua	109.686		5.0	13	100479	Type 48	*
	Drive CM Acqua	109.687		5.0	16	100479	Type 48	*
	Drive CM Acqua	109.898		5.0	18	100479	Type 48	*
	Drive CM Neoporos	109.1004		3.5	8	100479	Type 48	*
	Drive CM Neoporos	109.994		3.5	10	100479	Type 48	*
	Drive CM Neoporos	109.1005		3.5	11.5	100479	Type 48	*
	Drive CM Neoporos	109.995		3.5	13	100479	Type 48	*
	Drive CM Neoporos	109.996		3.5	16	100479	Type 48	*
	Drive CM Neoporos	109.1006		3.5	18	100479	Type 48	*
	Drive CM Neoporos	109.1001		4.3	8	100479	Type 48	*
	Drive CM Neoporos	109.1000		4.3	10	100479	Type 48	*
	Drive CM Neoporos	109.991		4.3	11.5	100479	Type 48	*
	Drive CM Neoporos	109.992		4.3	13	100479	Type 48	*
	Drive CM Neoporos	109.993		4.3	16	100479	Type 48	*
	Drive CM Neoporos	109.1007		4.3	18	100479	Type 48	*
	Drive CM Neoporos	109.1002		5.0	8	100479	Type 48	*
	Drive CM Neoporos	109.997		5.0	10	100479	Type 48	*
	Drive CM Neoporos	109.1003		5.0	11.5	100479	Type 48	*
	Drive CM Neoporos	109.998		5.0	13	100479	Type 48	*
	Drive CM Neoporos	109.999		5.0	16	100479	Type 48	*
	Drive CM Neoporos	109.1008		5.0	18	100479	Type 48	*
	Drive GM	109.958	GM	3.5	8	100543	Type 60	New
	Drive GM	109.959	GM	3.5	10	100543	Type 60	New
	Drive GM	109.960	GM	3.5	11.5	100543	Type 60	New
	Drive GM	109.961	GM	3.5	13	100543	Type 60	New
	Drive GM	109.962	GM	3.5	16	100543	Type 60	New
	Drive GM	109.963	GM	3.5	18	100543	Type 60	New
	Drive GM	109.964	GM	4.3	8	100543	Type 60	New
	Drive GM	109.965	GM	4.3	10	100543	Type 60	New
	Drive GM	109.966	GM	4.3	11.5	100543	Type 60	New
	Drive GM	109.967	GM	4.3	13	100543	Type 60	New
	Drive GM	109.968	GM	4.3	16	100543	Type 60	New
	Drive GM	109.969	GM	4.3	18	100543	Type 60	New



Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>NEODENT</b>	Drive GM	109.970	GM	5.0	8	100543	Type 60 New
	Drive GM	109.971	GM	5.0	10	100543	Type 60 New
	Drive GM	109.972	GM	5.0	11.5	100543	Type 60 New
	Drive GM	109.973	GM	5.0	13	100543	Type 60 New
	Drive GM	109.974	GM	5.0	16	100543	Type 60 New
	Drive GM	109.975	GM	5.0	18	100543	Type 60 New
	Drive Ti	109.757	4.3	4.3	8	100353	Type 01 *
	Drive Ti	109.758	4.3	4.3	10	100353	Type 01 *
	Drive Ti	109.759	4.3	4.3	11.5	100353	Type 01 *
	Drive Ti	109.760	4.3	4.3	13	100353	Type 01 *
	Drive Ti	109.761	4.3	4.3	16	100353	Type 01 *
	Drive Ti	109.762	5.0	5.0	8	100353	Type 01 *
	Drive Ti	109.763	5.0	5.0	10	100353	Type 01 *
	Drive Ti	109.764	5.0	5.0	11.5	100353	Type 01 *
	Drive Ti	109.765	5.0	5.0	13	100353	Type 01 *
	Drive Ti	109.766	5.0	5.0	16	100353	Type 01 *
	Helix GM	109.943	GM	3.5	8	100543	Type 60 New
	Helix GM	109.944	GM	3.5	10	100543	Type 60 New
	Helix GM	109.945	GM	3.5	11.5	100543	Type 60 New
	Helix GM	109.946	GM	3.5	13	100543	Type 60 New
	Helix GM	109.947	GM	3.5	16	100543	Type 60 New
	Helix GM	109.988	GM	3.5	18	100543	Type 60 New
	Helix GM	109.976	GM	3.75	8	100543	Type 60 New
	Helix GM	109.977	GM	3.75	10	100543	Type 60 New
	Helix GM	109.978	GM	3.75	11.5	100543	Type 60 New
	Helix GM	109.979	GM	3.75	13	100543	Type 60 New
	Helix GM	109.980	GM	3.75	16	100543	Type 60 New
	Helix GM	109.981	GM	3.75	18	100543	Type 60 New
	Helix GM	109.982	GM	4.0	8	100543	Type 60 New
	Helix GM	109.983	GM	4.0	10	100543	Type 60 New
	Helix GM	109.984	GM	4.0	11.5	100543	Type 60 New
	Helix GM	109.985	GM	4.0	13	100543	Type 60 New
	Helix GM	109.986	GM	4.0	16	100543	Type 60 New
	Helix GM	109.987	GM	4.0	18	100543	Type 60 New
	Helix GM	109.948	GM	4.3	8	100543	Type 60 New

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>NEODENT</b>	Helix GM	109.949	GM	4.3	10	100543	Type 60 New
	Helix GM	109.950	GM	4.3	11	100543	Type 60 New
	Helix GM	109.951	GM	4.3	13	100543	Type 60 New
	Helix GM	109.952	GM	4.3	16	100543	Type 60 New
	Helix GM	109.989	GM	4.3	18	100543	Type 60 New
	Helix GM	109.953	GM	5.0	8	100543	Type 60 New
	Helix GM	109.954	GM	5.0	10	100543	Type 60 New
	Helix GM	109.955	GM	5.0	11.5	100543	Type 60 New
	Helix GM	109.956	GM	5.0	13	100543	Type 60 New
	Helix GM	109.957	GM	5.0	16	100543	Type 60 New
	Helix GM	109.990	GM	5.0	18	100543	Type 60 New
	Helix GM	109.1009	GM	6.0	8	100543	Type 60 New
	Helix GM	109.1010	GM	6.0	10	100543	Type 60 New
	Helix GM	109.1011	GM	6.0	11.5	100543	Type 60 New
	Helix GM	109.1012	GM	6.0	13	100543	Type 60 New
	Titamax CM	109.613		3.5	7	100388	Type 16 *
	Titamax CM	109.614		3.5	8	100388	Type 16 *
	Titamax CM	109.615		3.5	9	100388	Type 16 *
	Titamax CM	109.616		3.5	11	100388	Type 16 *
	Titamax CM	109.617		3.5	13	100388	Type 16 *
	Titamax CM	109.618		3.5	15	100388	Type 16 *
	Titamax CM	109.619		3.5	17	100388	Type 16 *
	Titamax CM	109.606		3.75	7	100388	Type 16 *
	Titamax CM	109.607		3.75	8	100388	Type 16 *
	Titamax CM	109.608		3.75	9	100388	Type 16 *
	Titamax CM	109.609		3.75	11	100388	Type 16 *
	Titamax CM	109.610		3.75	13	100388	Type 16 *
	Titamax CM	109.611		3.75	15	100388	Type 16 *
	Titamax CM	109.612		3.75	17	100388	Type 16 *
	Titamax CM	109.630		4.0	7	100388	Type 16 *
	Titamax CM	109.631		4.0	8	100388	Type 16 *
	Titamax CM	109.632		4.0	9	100388	Type 16 *
	Titamax CM	109.633		4.0	11	100388	Type 16 *
	Titamax CM	109.620		4.0	13	100388	Type 16 *
	Titamax CM	109.634		4.0	15	100388	Type 16 *

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg	
<b>NEODENT</b>	Titamax CM	109.635		4.0	17	100388	Type 16	*
	Titamax CM	109.642		5.0	7	100388	Type 16	*
	Titamax CM	109.643		5.0	8	100388	Type 16	*
	Titamax CM	109.644		5.0	9	100388	Type 16	*
	Titamax CM	109.645		5.0	11	100388	Type 16	*
	Titamax CM	109.646		5.0	13	100388	Type 16	*
	Titamax CM EX	109.661		3.5	9	100388	Type 16	*
	Titamax CM EX	109.662		3.5	11	100388	Type 16	*
	Titamax CM EX	109.663		3.5	13	100388	Type 16	*
	Titamax CM EX	109.664		3.5	15	100388	Type 16	*
	Titamax CM EX	109.665		3.5	17	100388	Type 16	*
	Titamax CM EX	109.666		3.75	9	100388	Type 16	*
	Titamax CM EX	109.667		3.75	11	100388	Type 16	*
	Titamax CM EX	109.668		3.75	13	100388	Type 16	*
	Titamax CM EX	109.669		3.75	15	100388	Type 16	*
	Titamax CM EX	109.670		3.75	17	100388	Type 16	*
	Titamax CM EX	109.671		3.75	19	100388	Type 16	*
	Titamax CM EX	109.636		4.0	9	100388	Type 16	*
	Titamax CM EX	109.637		4.0	11	100388	Type 16	*
	Titamax CM EX	109.638		4.0	13	100388	Type 16	*
	Titamax CM EX	109.639		4.0	15	100388	Type 16	*
	Titamax CM EX	109.640		4.0	17	100388	Type 16	*
	Titamax CM EX	109.641		4.0	19	100388	Type 16	*
	Titamax GM	109.906	GM	3.5	7	100543	Type 60	New
	Titamax GM	109.907	GM	3.5	8	100543	Type 60	New
	Titamax GM	109.908	GM	3.5	9	100543	Type 60	New
	Titamax GM	109.909	GM	3.5	11	100543	Type 60	New
	Titamax GM	109.910	GM	3.5	13	100543	Type 60	New
	Titamax GM	109.911	GM	3.5	15	100543	Type 60	New
	Titamax GM	109.912	GM	3.5	17	100543	Type 60	New
	Titamax GM	109.899	GM	3.75	7	100543	Type 60	New
	Titamax GM	109.900	GM	3.75	8	100543	Type 60	New
	Titamax GM	109.901	GM	3.75	9	100543	Type 60	New
	Titamax GM	109.902	GM	3.75	11	100543	Type 60	New
	Titamax GM	109.903	GM	3.75	13	100543	Type 60	New

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg	
<b>NEODENT</b>	Titamax GM	109.904	GM	3.75	15	100543	Type 60	New
	Titamax GM	109.905	GM	3.75	17	100543	Type 60	New
	Titamax GM	109.913	GM	4.0	7	100543	Type 60	New
	Titamax GM	109.914	GM	4.0	8	100543	Type 60	New
	Titamax GM	109.915	GM	4.0	9	100543	Type 60	New
	Titamax GM	109.916	GM	4.0	11	100543	Type 60	New
	Titamax GM	109.917	GM	4.0	13	100543	Type 60	New
	Titamax GM	109.918	GM	4.0	15	100543	Type 60	New
	Titamax GM	109.919	GM	4.0	17	100543	Type 60	New
	Titamax GM	109.920	GM	5.0	7	100543	Type 60	New
	Titamax GM	109.921	GM	5.0	8	100543	Type 60	New
	Titamax GM	109.922	GM	5.0	9	100543	Type 60	New
	Titamax GM	109.923	GM	5.0	11	100543	Type 60	New
	Titamax GM	109.924	GM	5.0	13	100543	Type 60	New
	Titamax GT			3.7		100368	Type 09	
	Titamax GT			4.0		100368	Type 09	
	Titamax GT			4.5		100368	Type 09	
	Titamax Ti	109.284	4.1	3.75	9	100353	Type 01	*
	Titamax Ti	109.285	4.1	3.75	11	100353	Type 01	*
	Titamax Ti	109.286	4.1	3.75	13	100353	Type 01	*
	Titamax Ti	109.287	4.1	3.75	15	100353	Type 01	*
	Titamax Ti	109.288	4.1	3.75	17	100353	Type 01	*
	Titamax Ti	109.302	4.1	4.0	9	100353	Type 01	*
	Titamax Ti	109.303	4.1	4.0	11	100353	Type 01	*
	Titamax Ti	109.304	4.1	4.0	13	100353	Type 01	*
	Titamax Ti	109.305	4.1	4.0	15	100353	Type 01	*
	Titamax Ti	109.306	4.1	4.0	17	100353	Type 01	*
	Titamax Ti	109.327	5.0	4.5	9	100357	Type 03	*
	Titamax Ti	109.328	5.0	4.5	11	100357	Type 03	*
	Titamax Ti	109.329	5.0	4.5	13	100357	Type 03	*
	Titamax Ti	109.294	5.0	5.0	7	100357	Type 03	*
	Titamax Ti	109.295	5.0	5.0	9	100357	Type 03	*
	Titamax Ti	109.296	5.0	5.0	11	100357	Type 03	*
	Titamax Ti	109.297	5.0	5.0	13	100357	Type 03	*
	Titamax Ti EX	109.582	4.1	3.75	9	100353	Type 01	*

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg	
<b>NEODENT</b>	Titamax Ti EX	109.583	4.1	3.75	11	100353	Type 01	*
	Titamax Ti EX	109.584	4.1	3.75	13	100353	Type 01	*
	Titamax Ti EX	109.585	4.1	3.75	15	100353	Type 01	*
	Titamax Ti EX	109.586	4.1	3.75	17	100353	Type 01	*
	Titamax Ti EX	109.587	4.1	3.75	19	100353	Type 01	*
	Titamax Ti EX	109.588	4.1	4.0	9	100353	Type 01	*
	Titamax Ti EX	109.589	4.1	4.0	11	100353	Type 01	*
	Titamax Ti EX	109.590	4.1	4.0	13	100353	Type 01	*
	Titamax Ti EX	109.591	4.1	4.0	15	100353	Type 01	*
	Titamax Ti EX	109.592	4.1	4.0	17	100353	Type 01	*
	Titamax Ti EX	109.593	4.1	4.0	19	100353	Type 01	*
	Titamax II Plus			3.75		100402	Type 21	
	Titamax II Plus			5.0		100402	Type 21	
<b>NEOSS</b>	Neoss Implant			3.25		100455	Type 38	
	ProActive Straight			3.5		100398	Type 18	
	ProActive Straight			4.0		100398	Type 18	
	ProActive Straight			4.5		100398	Type 18	
	ProActive Straight			5.0		100398	Type 18	
	ProActive Straight			5.5		100398	Type 18	
	ProActive Straight			6.0		100398	Type 18	
	ProActive Tapered			3.5		100398	Type 18	
	ProActive Tapered			4.0		100398	Type 18	
	ProActive Tapered			4.5		100398	Type 18	
	ProActive Tapered			5.0		100398	Type 18	
	ProActive Tapered			5.5		100398	Type 18	
	ProActive Tapered			6.0		100398	Type 18	
<b>NIPPON PISTON RING</b>	1 Stage Screw RH			3.3		100398	Type 18	
	1 Stage Screw RH			4.0		100398	Type 18	
	1 Stage Screw RL			3.3		100398	Type 18	
	1 Stage Screw RL			4.0		100398	Type 18	
	1 Stage Screw WH			5.0		100398	Type 18	
	1 Stage Screw WL			5.0		100398	Type 18	
	2 Stage Cylinder			3.3		100398	Type 18	

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>NIPPON PISTON RING</b>	2 Stage Cylinder			4.0		100398	Type 18
	2 Stage Cylinder			5.0		100398	Type 18
	2 Stage Screw			3.3		100398	Type 18
	2 Stage Screw			4.0		100398	Type 18
	2 Stage Screw			5.0		100398	Type 18
<b>NOBEL BIOCARE</b>	Brånemark system		NP			100355	Type 02
	Brånemark system		RP			100353	Type 01
	Brånemark system		WP			100357	Type 03
	Brånemark system Novum					100353	Type 01
	Brånemark system Shorty		RP			100353	Type 01
	Brånemark system Shorty		WP			100357	Type 03
	NobelActive		NP	3.5		100543	Type 60
	NobelActive		RP	4.3		100544	Type 61
	NobelActive		RP	5.0		100544	Type 61
	NobelActive		WP	5.5		100350	Type 04
	NobelActive			3.0		100545	Type 62
	NobelDirect		NP			100366	Type 08
	NobelDirect		RP			100368	Type 09
	NobelDirect		WP			100368	Type 09
	NobelParallel Conical Connection		NP	3.75		100543	Type 60
	NobelParallel Conical Connection		RP	4.3		100544	Type 61
	NobelParallel Conical Connection		RP	5.0		100544	Type 61
	NobelParallel Conical Connection		WP	5.5		100350	Type 04
	NobelReplace Conical Connection		NP	3.5		100543	Type 60
	NobelReplace Conical Connection		RP	4.3		100544	Type 61
	NobelReplace Conical Connection		RP	5.0		100544	Type 61
	NobelReplace Conical Connection PMC		NP			100543	Type 60
	NobelReplace Conical Connection PMC		RP			100544	Type 61
	NobelReplace Plateforme Shift		NP			100359	Type 12
	NobelReplace Plateforme Shift		RP			100361	Type 13
	NobelReplace Plateforme Shift		WP			100361	Type 13
	NobelReplace Straight Groovy		NP			100359	Type 12
	NobelReplace Straight Groovy		RP			100361	Type 13
	NobelReplace Straight Groovy		WP			100361	Type 13

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>NOBEL BIOCARE</b>	NobelReplace Tapered Groovy		NP			100359	Type 12
	NobelReplace Tapered Groovy		RP			100361	Type 13
	NobelReplace Tapered Groovy		WP			100361	Type 13
	NobelReplace Tapered Groovy			6.0		100462	Type 40
	NobelSpeedy Groovy		NP			100355	Type 02
	NobelSpeedy Groovy		RP			100353	Type 01
	NobelSpeedy Groovy		WP			100357	Type 03
	NobelSpeedy Replace		NP			100359	Type 12
	NobelSpeedy Replace		RP			100361	Type 13
	NobelSpeedy Replace		WP			100361	Type 13
	NobelSpeedy Replace			6.0		100462	Type 40
	NobelSpeedy Shorty		RP			100353	Type 01
	NobelSpeedy Shorty		WP			100357	Type 03
	Replace Select Straight		NP			100359	Type 12
	Replace Select Straight		RP			100361	Type 13
	Replace Select Straight		WP			100361	Type 13
	Replace Select Tapered		NP			100359	Type 12
	Replace Select Tapered		RP			100361	Type 13
	Replace Select Tapered		WP			100361	Type 13
	Replace Select Tapered			6.0		100462	Type 40
	Replace Select Tapered PMC		NP			100359	Type 12
	Replace Select Tapered PMC		RP			100361	Type 13
	Replace Select Tapered PMC		WP			100361	Type 13
	Replace Select Tapered PMC			6.0		100462	Type 40
	Replace Select TC		NP			100359	Type 12
	Replace Select TC		RP			100361	Type 13
<b>NUCLEOSS</b>	Nucleoss T3			3.4		100376	Type 05
	Nucleoss T3			3.8		100376	Type 05
	Nucleoss T3			4.2		100376	Type 05
	Nucleoss T3			5.0		100376	Type 05
	Nucleoss T4			3.4		100376	Type 05
	Nucleoss T4			3.8		100376	Type 05
	Nucleoss T4			4.2		100376	Type 05
	Nucleoss T5			5.0		100380	Type 07

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>OCO BIOMEDICAL</b>	Engage		3.0		3.2	100455	Type 38
	Engage		3.5		4.0	100425	Type 26
	Engage		4.5		5.0	100425	Type 26
	Engage		4.5		6.0	100431	Type 27
	ERI				4.0	100380	Type 07
	ERI				5.0	100380	Type 07
	Macro				6.0	100478	Type 47
	TSI				3.25	100446	Type 35
	TSI				4.0	100380	Type 07
	TSI				5.0	100380	Type 07
<b>OSSEOLINK</b>					3.5	100402	Type 21
					4.0	100402	Type 21
					5.0	100402	Type 21
<b>OSSTEM</b>	GSII		Mini			100376	Type 05
	GSII		Regular			100378	Type 06
	GSII		Wide			100378	Type 06
	GSII		Ultra Wide			100378	Type 06
	GSIII		Mini			100376	Type 05
	GSIII		Regular			100378	Type 06
	SSI		Regular			100350	Type 04
	SSI		Wide			100350	Type 04
	SSII		Regular			100350	Type 04
	SSII		Wide			100350	Type 04
	SSIII		Regular			100350	Type 04
	SSIII		Wide			100350	Type 04
	TSIII CA				3.5	100376	Type 05
	TSIII SA				3.0	100372	Type 11
	TSIII SA				4.0	100478	Type 47
	TSIII SA				4.5	100478	Type 47
	TSIII SA				5.0	100478	Type 47
	TSIII SA Fixture – NoMount				4.5	100478	Type 47
	TSIII SA Ultra-Wide				6.0	100350	Type 04
	TSIII SA Ultra-Wide				7.0	100350	Type 04



Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>OSSTEM</b>	USII		Mini			100355	Type 02
	USII		Regular			100353	Type 01
	USII		Wide			100357	Type 03
	USII		Ultra Wide			100357	Type 03
	USIII		Mini			100355	Type 02
	USIII		Regular			100353	Type 01
	USIII		Wide			100357	Type 03
	USIII		Ultra Wide			100357	Type 03
	USIV Conical		Mini			100366	Type 08
<b>OSTEOPLANT</b>	Two-piece					100370	Type 10
<b>OSTEO-TI</b>	Combi Dental Implant			3.25		100425	Type 26
	Combi Dental Implant			3.75		100425	Type 26
	Combi Dental Implant			4.5		100425	Type 26
	Combi Dental Implant			5.5		100425	Type 26
<b>OTICON MEDICAL</b>	Ponto Implant			3.0		100368	Type 09
	Ponto Implant			4.0		100368	Type 09
	Ponto Implant		Wide			100368	Type 09
<b>PALTOP</b>	PALTOP		NP			100455	Type 38
	PALTOP		SP			100480	Type 49
<b>PHIBO</b>	Aurea EVO		NP			100554	Type 63
	Aurea EVO		RP			100425	Type 26
	Aurea EVO		WP			100425	Type 26
	Aurea		NP			100554	Type 63
	Aurea		RP			100431	Type 27
	Aurea		WP			100440	Type 32
	BNT				S3	100353	Type 01
	BNT				S4	100353	Type 01
	BNT				S5	100353	Type 01
	TSA				S3	100455	Type 38
	TSA				S4	100431	Type 27

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>PHIBO</b>	TSA				S5	100431	Type 27
	TSA Advance				S3	100455	Type 38
	TSA Advance				S4	100431	Type 27
	TSA Advance				S5	100431	Type 27
	TSH				S3	100353	Type 01
	TSH				S4	100353	Type 01
	TSH				S5	100353	Type 01
<b>P-I BRÄNEMARK PHILOSOPHY</b>	Internal Grip Amplified				3.5	100554	Type 63
	Internal Grip Amplified				4.3	100482	Type 51
	Internal Grip Amplified				5.1	100482	Type 51
	Internal Grip External Hex				3.5	100455	Type 38
	Internal Grip External Hex				3.75	100398	Type 18
	Internal Grip External Hex				4.1	100398	Type 18
	Morse Taper				3.3	100402	Type 21
	Morse Taper				3.75	100402	Type 21
	Morse Taper				4.8	100402	Type 21
<b>PLAN 1 HEALTH</b>	GTB				3.3	100497	Type 57
	GTB				3.6	100497	Type 57
	GTB				4.3	100497	Type 57
	Hexcel EHC, External Hex					100353	Type 01
	Hexcel IHC, Internal Hex					100425	Type 26
	Hexcel IHT, Internal Hex					100425	Type 26
<b>PLATON</b>	EL Implant Narrow				3.7	100543	Type 60
	EL Implant Regular				5.1	100478	Type 47
	Type I					100350	Type 04
	Type II					100350	Type 04
	Type IV Pro					100478	Type 47
<b>PRODENT ITALIA</b>	Perform SM				3.6	100480	Type 49
	Perform SM				4.0	100480	Type 49
	Perform SM				4.5	100380	Type 07
	Perform SM				5.0	100380	Type 07

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg	
<b>PW PLUS</b>	Dif				3.8	100480	Type 49	
	Dif				4.2	100380	Type 07	
<b>REVIOS</b>					3.8	100425	Type 26	
					4.3	100425	Type 26	
					5.0	100425	Type 26	
<b>REX IMPLANTS</b>	Rex TL				1.8	100350	Type 04	
	Rex TL				2.9	100480	Type 49	
<b>RUETSCHI</b>	BL				4.1	100670	Type 80	New
	TLR				4.7	100670	Type 80	New
	TLW				4.7	100670	Type 80	New
<b>S.I.N IMPLANT SYSTEM</b>	External Hex				3.4	100353	Type 01	
	External Hex				4.1	100353	Type 01	
	External Hex				5.0	100353	Type 01	
	Internal Hex				3.8	100386	Type 15	
	Internal Hex				4.5	100386	Type 15	
	Internal Hex				5.5	100386	Type 15	
	Strong SW HI				4.5	100479	Type 48	
	Strong SW Morse				4.5	100484	Type 53	
<b>SDS</b>	SDS 2.0				3.8	100589	Type 69	
<b>SERVO-DENTAL</b>	Trias				3.8	100378	Type 06	
<b>SHINHUNG</b>	Luna		NP		3.5	100376	Type 05	
	Luna		RP		4.0	100478	Type 47	
	Luna		RP		4.5	100478	Type 47	
	Luna		RP		5.0	100478	Type 47	
	Sola		NP		3.5	100355	Type 02	
	Sola		RP		4.0	100353	Type 01	
	Sola		RP		4.5	100353	Type 01	
	Sola		WP		5.0	100357	Type 03	

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>SHINHUNG</b>	Stella		RP		4.0	100544	Type 61
	Stella		RP		4.5	100544	Type 61
	Stella		WP		4.5	100544	Type 61
	Stella		WP		5.0	100544	Type 61
<b>SHOFU</b>	BFIMP30-OSTL				3.0	100545	Type 62
	BFIMP35-OSTL				3.5	100404	Type 22
	BFIMP40-OSTL				4.0	100543	Type 60
<b>SIC IMPLANT SYSTEM</b>	Ace				3.4	100466	Type 44
	Ace				4.0	100466	Type 44
	Ace				5.0	100466	Type 44
<b>SOUTHERN IMPLANTS</b>	External Hex					100353	Type 01
	Extra-oral			3.75		100353	Type 01
	IT Connection					100350	Type 04
	MAX				7.0	100353	Type 01
	MAX				8.0	100353	Type 01
	MAX				9.0	100353	Type 01
	MAX-IT				7.0	100380	Type 07
	MAX-IT				8.0	100380	Type 07
	MAX-IT				9.0	100380	Type 07
	Tri-MAX				7.0	100412	Type 24
	Tri-MAX				8.0	100462	Type 40
	Tri-MAX				9.0	100462	Type 40
	Tri-Nex				3.5	100359	Type 12
	Tri-Nex				4.3	100361	Type 13
Tri-Nex				5.0	100361	Type 13	
<b>STERNGOLD DENTAL</b>	External Hexagon		RP		3.75	100353	Type 01
	External Hexagon		RP		4.0	100353	Type 01
	External Hexagon		RP		5.0	100353	Type 01
	Stern IC Internal Octagon		RN		4.1	100350	Type 04
<b>STRAUMANN</b>	Bone Level Implants Roxolid SLActive		NC		3.3	100484	Type 53

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>STRAUMANN</b>	Bone Level Implants Roxolid SLActive		RC		4.1	100485	Type 54
	Bone Level Implants Roxolid SLActive		RC		4.8	100485	Type 54
	Bone Level Implants Titanium SLA		NC		3.3	100484	Type 53
	Bone Level Implants Titanium SLA		RC		4.1	100485	Type 54
	Bone Level Implants Titanium SLA		RC		4.8	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.3508	NC		3.3 8	100484	Type 53
	Bone Level Tapered Implants Roxolid SLA	021.3510	NC		3.3 10	100484	Type 53
	Bone Level Tapered Implants Roxolid SLA	021.3512	NC		3.3 12	100484	Type 53
	Bone Level Tapered Implants Roxolid SLA	021.3514	NC		3.3 14	100484	Type 53
	Bone Level Tapered Implants Roxolid SLA	021.3516	NC		3.3 16	100484	Type 53
	Bone Level Tapered Implants Roxolid SLA	021.3518	NC		3.3 18	100484	Type 53
	Bone Level Tapered Implants Roxolid SLA	021.5508	RC		4.1 8	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.5510	RC		4.1 10	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.5512	RC		4.1 12	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.5514	RC		4.1 14	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.5516	RC		4.1 16	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.5518	RC		4.1 18	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.7508	RC		4.8 8	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.7510	RC		4.8 10	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.7512	RC		4.8 12	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.7514	RC		4.8 14	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.7516	RC		4.8 16	100485	Type 54
	Bone Level Tapered Implants Roxolid SLA	021.7518	RC		4.8 18	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.3308	NC		3.3 8	100484	Type 53
	Bone Level Tapered Implants Roxolid SLActive	021.3310	NC		3.3 10	100484	Type 53
	Bone Level Tapered Implants Roxolid SLActive	021.3312	NC		3.3 12	100484	Type 53
	Bone Level Tapered Implants Roxolid SLActive	021.3314	NC		3.3 14	100484	Type 53
	Bone Level Tapered Implants Roxolid SLActive	021.3316	NC		3.3 16	100484	Type 53
	Bone Level Tapered Implants Roxolid SLActive	021.3318	NC		3.3 18	100484	Type 53
	Bone Level Tapered Implants Roxolid SLActive	021.5308	RC		4.1 8	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.5310	RC		4.1 10	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.5312	RC		4.1 12	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.5314	RC		4.1 14	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.5316	RC		4.1 16	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.5318	RC		4.1 18	100485	Type 54

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>STRAUMANN</b>	Bone Level Tapered Implants Roxolid SLActive	021.7308	RC	4.8	8	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.7310	RC	4.8	10	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.7312	RC	4.8	12	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.7314	RC	4.8	14	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.7316	RC	4.8	16	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive	021.7318	RC	4.8	18	100485	Type 54
	Bone Level Tapered Implants Roxolid SLActive		SC	2.9		100404	Type 22
	Bone Level Tapered Implants Titanium SLA	021.3408	NC	3.3	8	100484	Type 53
	Bone Level Tapered Implants Titanium SLA	021.3410	NC	3.3	10	100484	Type 53
	Bone Level Tapered Implants Titanium SLA	021.3412	NC	3.3	12	100484	Type 53
	Bone Level Tapered Implants Titanium SLA	021.3414	NC	3.3	14	100484	Type 53
	Bone Level Tapered Implants Titanium SLA	021.3416	NC	3.3	16	100484	Type 53
	Bone Level Tapered Implants Titanium SLA	021.3418	NC	3.3	18	100484	Type 53
	Bone Level Tapered Implants Titanium SLA	021.5408	RC	4.1	8	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.5410	RC	4.1	10	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.5412	RC	4.1	12	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.5414	RC	4.1	14	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.5416	RC	4.1	16	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.5418	RC	4.1	18	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.7408	RC	4.8	8	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.7410	RC	4.8	10	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.7412	RC	4.8	12	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.7414	RC	4.8	14	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.7416	RC	4.8	16	100485	Type 54
	Bone Level Tapered Implants Titanium SLA	021.7418	RC	4.8	18	100485	Type 54
	Extra Oral					100353	Type 01
	Narrow Neck					100390	Type 17
	Ortho					100368	Type 09
	Palatal					100368	Type 09
	PURE Ceramic ZLA, ZrO2	032.000S	RD	4.1	8	100671	Type 81 New
	PURE Ceramic ZLA, ZrO2	032.001S	RD	4.1	10	100671	Type 81 New
	PURE Ceramic ZLA, ZrO2	032.002S	RD	4.1	12	100671	Type 81 New
	PURE Ceramic ZLA, ZrO2	032.003S	RD	4.1	14	100671	Type 81 New
	Standard Implants Roxolid SLActive		RN	3.3		100350	Type 04
	Standard Implants Roxolid SLActive		RN	4.1		100350	Type 04

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>STRAUMANN</b>	Standard Implants Roxolid SLActive		RN	4.8		100350	Type 04
	Standard Implants Roxolid SLActive		WN	4.8		100350	Type 04
	Standard Implants Titanium SLA		RN	3.3		100350	Type 04
	Standard Implants Titanium SLA		RN	4.1		100350	Type 04
	Standard Implants Titanium SLA		RN	4.8		100350	Type 04
	Standard Implants Titanium SLA		WN	4.8		100350	Type 04
	Standard Plus Implants Roxolid SLActive		NNC	3.3		100560	Type 68
	Standard Plus Implants Roxolid SLActive		RN	3.3		100350	Type 04
	Standard Plus Implants Roxolid SLActive		RN	4.1		100350	Type 04
	Standard Plus Implants Roxolid SLActive		RN	4.8		100350	Type 04
	Standard Plus Implants Roxolid SLActive		WN	4.8		100350	Type 04
	Standard Plus Implants Titanium SLA		RN	3.3		100350	Type 04
	Standard Plus Implants Titanium SLA		RN	4.1		100350	Type 04
	Standard Plus Implants Titanium SLA		RN	4.8		100350	Type 04
	Standard Plus Implants Titanium SLA		WN	4.8		100350	Type 04
	Tapered Effect Implants Roxolid SLActive		RN	3.3		100350	Type 04
	Tapered Effect Implants Roxolid SLActive		RN	4.1		100350	Type 04
	Tapered Effect Implants Roxolid SLActive		WN	4.8		100350	Type 04
	Tapered Effect Implants Titanium SLA		RN	3.3		100350	Type 04
	Tapered Effect Implants Titanium SLA		RN	4.1		100350	Type 04
	Tapered Effect Implants Titanium SLA		WN	4.8		100350	Type 04
<b>SWEDEN &amp; MARTINA</b>	EZ					100380	Type 07
	Global			3.8		100479	Type 48
	Global			4.3		100479	Type 48
	Global			4.8		100479	Type 48
	Global			5.5		100479	Type 48
	Kohno			3.8		100440	Type 32
	Kohno			4.25		100380	Type 07
	Kohno			5.0		100380	Type 07
	Premium - new			3.3		100440	Type 32
	Premium - new			3.8		100440	Type 32
	Premium			4.25		100380	Type 07
	Premium			5.0		100380	Type 07
	Premium - old			3.3		100402	Type 21

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>SYBRON</b>	Endopore External				3.5	100436	Type 30
	Endopore External				4.1	100353	Type 01
	Endopore External				5.0	100353	Type 01
	Entegra External				3.25	100353	Type 01
	Entegra External				4.0	100353	Type 01
	Pitt-Easy				4.0	100412	Type 24
	Pitt-Easy				4.9	100412	Type 24
	SybronPRO XRT				3.3	100378	Type 06
	SybronPRO XRT				3.5	100378	Type 06
	SybronPRO XRT				4.1	100378	Type 06
<b>TA-DENT</b>					3.3	100455	Type 38
					4.8	100378	Type 06
					5.6	100480	Type 49
<b>TAG DENTAL</b>	Axis				3.3	100431	Type 27
	Axis				3.75	100431	Type 27
	Axis				4.2	100431	Type 27
	Axis				5.0	100431	Type 27
	Axis				6.0	100431	Type 27
	Massif				3.75	100431	Type 27
	Massif				4.2	100431	Type 27
	Massif				5.0	100431	Type 27
	Massif				6.0	100431	Type 27
<b>TBR</b>	External Hexagon					100353	Type 01
	Internal Octagon					100398	Type 18
	M-implant				3.2	100402	Type 21
	M-implant				3.9	100402	Type 21
	M-implant				4.7	100402	Type 21
	Z1-M implant				3.2	100402	Type 21
	Z1-M implant				3.9	100402	Type 21
	Z1-M implant				4.7	100402	Type 21
<b>TEKKA</b>	Inkone				3.5	100480	Type 49



Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>TEKKA</b>	Inkone				4.5	100480	Type 49
	Krestal				3.5	100378	Type 06
	Krestal				4.0	100378	Type 06
	Krestal				4.2	100378	Type 06
	Krestal				5.1	100378	Type 06
	Progress				3.4	100378	Type 06
	Progress				4.0	100378	Type 06
	T-quest				4.0	100480	Type 49
	T-quest				4.5	100480	Type 49
	Twinkone				3.5	100390	Type 17
<b>TFI SYSTEM</b>	TFI				3.3	100425	Type 26
	TFI				3.75	100425	Type 26
	TFI				4.2	100425	Type 26
	TFI				5.0	100425	Type 26
<b>THOMMEN MEDICAL</b>	SPI System Plateforme				3.5	100404	Type 22
	SPI System Plateforme				4.0	100402	Type 21
	SPI System Plateforme				4.5	100386	Type 15
	SPI System Plateforme				5.0	100386	Type 15
	SPI System Plateforme				6.0	100386	Type 15
<b>TIDAL SPIRAL</b>					3.75	100431	Type 27
					4.2	100431	Type 27
					5.0	100431	Type 27
					6.0	100431	Type 27
<b>TRATE</b>	Rootform				4.0	100497	Type 57
<b>TRI DENTAL IMPLANTS</b>	TRI-Octa 4.8				3.75	100544	Type 61
	TRI-Octa 4.8				4.1	100544	Type 61
	TRI-Octa 4.8				4.7	100544	Type 61
	TRI-Vent				3.5	100440	Type 32
<b>U-IMPL</b>	External Hex					100353	Type 01

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg	
<b>U-IMPL</b>	Internal Hex					100440	Type 32	
<b>WARANTEC</b>	Oneplant Implant			4.3		100412	Type 24	
	Oneplant Internal Tapered					100378	Type 06	
<b>WEGO</b>				3.4	9	100484	Type 53	
				4.3	7	100455	Type 38	
				4.3	9	100411	Type 23	
<b>WIELAND IMPLANT SYSTEM</b>				3.5		100472	Type 45	
				4.3		100472	Type 45	
				5.0		100472	Type 45	
<b>YAMAHACHI DENTAL</b>	μ-one HA Implant 2 piece			3.4		100455	Type 38	
	μ-one HA Implant 2 piece			3.8		100455	Type 38	
	μ-one HA Implant 2 piece			4.3		100431	Type 27	
	μ-one HA Implant 2 piece			5.0		100398	Type 18	
<b>ZERAMEX</b>	P6		RN			100674	Type 82	New
<b>ZIACOM</b>	Implante Zinic	ZSS3785	RP	3.7	8.5	100440	Type 32	New
	Implante Zinic	ZSS3710	RP	3.7	10	100440	Type 32	New
	Implante Zinic	ZSS3711	RP	3.7	11.5	100440	Type 32	New
	Implante Zinic	ZSS3713	RP	3.7	13	100440	Type 32	New
	Implante Zinic	ZSS3714	RP	3.7	14.5	100440	Type 32	New
	Implante Zinic	ZSS4085	RP	4.0	8.5	100440	Type 32	New
	Implante Zinic	ZSS4010	RP	4.0	10	100440	Type 32	New
	Implante Zinic	ZSS4011	RP	4.0	11.5	100440	Type 32	New
	Implante Zinic	ZSS4013	RP	4.0	13	100440	Type 32	New
	Implante Zinic	ZSS4014	RP	4.0	14.5	100440	Type 32	New
	Implante Zinic	ZSS4385	RP	4.3	8.5	100440	Type 32	New
	Implante Zinic	ZSS4310	RP	4.3	10	100440	Type 32	New
	Implante Zinic	ZSS4311	RP	4.3	11.5	100440	Type 32	New
	Implante Zinic	ZSS4313	RP	4.3	13	100440	Type 32	New
	Implante Zinic	ZSS4314	RP	4.3	14.5	100440	Type 32	New

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>ZIACOM</b>	Implante Zinic MT	ZSS3685M	RP	3.6	8.5	100440	Type 32 New
	Implante Zinic MT	ZSS3610M	RP	3.6	10	100440	Type 32 New
	Implante Zinic MT	ZSS3611M	RP	3.6	11.5	100440	Type 32 New
	Implante Zinic MT	ZSS3613M	RP	3.6	13	100440	Type 32 New
	Implante Zinic MT	ZSS4085M	RP	4.0	8.5	100440	Type 32 New
	Implante Zinic MT	ZSS4010M	RP	4.0	10	100440	Type 32 New
	Implante Zinic MT	ZSS4011M	RP	4.0	11.5	100440	Type 32 New
	Implante Zinic MT	ZSS4013M	RP	4.0	13	100440	Type 32 New
	Implante Zinic MT	ZSS4485M	RP	4.4	8.5	100440	Type 32 New
	Implante Zinic MT	ZSS4410M	RP	4.4	10	100440	Type 32 New
	Implante Zinic MT	ZSS4411M	RP	4.4	11.5	100440	Type 32 New
	Implante Zinic MT	ZSS4413M	RP	4.4	13	100440	Type 32 New
	Implante Zinic Shorty	ZSS4756HT	RP	4.75	6	100440	Type 32 New
	Implante Zinic Shorty	ZSS4757HT	RP	4.75	7	100440	Type 32 New
	Implante ZM1	ZM13685	RP	3.6	8.5	100353	Type 1 New
	Implante ZM1	ZM13610	RP	3.6	10	100353	Type 1 New
	Implante ZM1	ZM13611	RP	3.6	11.5	100353	Type 1 New
	Implante ZM1	ZM13613	RP	3.6	13	100353	Type 1 New
	Implante ZM1	ZM14085	RP	4.0	8.5	100353	Type 1 New
	Implante ZM1	ZM14010	RP	4.0	10	100353	Type 1 New
	Implante ZM1	ZM14011	RP	4.0	11.5	100353	Type 1 New
	Implante ZM1	ZM14013	RP	4.0	13	100353	Type 1 New
	Implante ZM1	ZM14485	RP	4.4	8.5	100353	Type 1 New
	Implante ZM1	ZM14410	RP	4.4	10	100353	Type 1 New
	Implante ZM1	ZM14411	RP	4.4	11.5	100353	Type 1 New
	Implante ZM1	ZM14413	RP	4.4	13	100353	Type 1 New
	Implante ZM4	ZM43785	RP	3.7	8.5	100353	Type 1 New
	Implante ZM4	ZM43710	RP	3.7	10	100353	Type 1 New
	Implante ZM4	ZM43711	RP	3.7	11.5	100353	Type 1 New
	Implante ZM4	ZM43713	RP	3.7	13	100353	Type 1 New
	Implante ZM4	ZM43714	RP	3.7	14.5	100353	Type 1 New
	Implante ZM4	ZM44085	RP	4.0	8.5	100353	Type 1 New
	Implante ZM4	ZM44010	RP	4.0	10	100353	Type 1 New
	Implante ZM4	ZM44011	RP	4.0	11.5	100353	Type 1 New
	Implante ZM4	ZM44013	RP	4.0	13	100353	Type 1 New

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>ZIACOM</b>	Implante ZM4	ZM44014	RP		4.0 14.5	100353	Type 1 New
	Implante ZM4	ZM44385	RP		4.3 8.5	100353	Type 1 New
	Implante ZM4	ZM44310	RP		4.3 10	100353	Type 1 New
	Implante ZM4	ZM44311	RP		4.3 11.5	100353	Type 1 New
	Implante ZM4	ZM44313	RP		4.3 13	100353	Type 1 New
	Implante ZM4	ZM44314	RP		4.3 14.5	100353	Type 1 New
	Implante ZM8	ZM83709N	RP		3.75 9	100353	Type 1 New
	Implante ZM8	ZM83711N	RP		3.75 11	100353	Type 1 New
	Implante ZM8	ZM83713N	RP		3.75 13	100353	Type 1 New
	Implante ZM8	ZM83715N	RP		3.75 15	100353	Type 1 New
	Implante ZM8	ZM84009N	RP		4.0 9	100353	Type 1 New
	Implante ZM8	ZM84011N	RP		4.0 11	100353	Type 1 New
	Implante ZM8	ZM84013N	RP		4.0 13	100353	Type 1 New
	Implante ZM8	ZM84015N	RP		4.0 15	100353	Type 1 New
	Implante ZM8	ZM84309N	RP		4.3 9	100353	Type 1 New
	Implante ZM8	ZM84311N	RP		4.3 11	100353	Type 1 New
	Implante ZM8	ZM84313N	RP		4.3 13	100353	Type 1 New
	Implante ZM8	ZM84315N	RP		4.3 15	100353	Type 1 New
	Implante ZM8	ZM83709S	RP		3.75 9	100353	Type 1 New
	Implante ZM8	ZM83711S	RP		3.75 11	100353	Type 1 New
	Implante ZM8	ZM83713S	RP		3.75 13	100353	Type 1 New
	Implante ZM8	ZM83715S	RP		3.75 15	100353	Type 1 New
	Implante ZM8	ZM84009S	RP		4.0 9	100353	Type 1 New
	Implante ZM8	ZM84011S	RP		4.0 11	100353	Type 1 New
	Implante ZM8	ZM84013S	RP		4.0 13	100353	Type 1 New
	Implante ZM8	ZM84015S	RP		4.0 15	100353	Type 1 New
	Implante ZM8	ZM84309S	RP		4.3 9	100353	Type 1 New
	Implante ZM8	ZM84311S	RP		4.3 11	100353	Type 1 New
	Implante ZM8	ZM84313S	RP		4.3 13	100353	Type 1 New
	Implante ZM8	ZM84315S	RP		4.3 15	100353	Type 1 New
	Implante ZM4 MT	ZM43685M	RP		3.6 8.5	100353	Type 1 New
	Implante ZM4 MT	ZM43610M	RP		3.6 10	100353	Type 1 New
	Implante ZM4 MT	ZM43611M	RP		3.6 11.5	100353	Type 1 New
	Implante ZM4 MT	ZM43613M	RP		3.6 13	100353	Type 1 New
	Implante ZM4 MT	ZM44085M	RP		4.0 8.5	100353	Type 1 New

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>ZIACOM</b>	Implante ZM4 MT	ZM44010M	RP	4.0	10	100353	Type 1 New
	Implante ZM4 MT	ZM44011M	RP	4.0	11.5	100353	Type 1 New
	Implante ZM4 MT	ZM44013M	RP	4.0	13	100353	Type 1 New
	Implante ZM4 MT	ZM44485M	RP	4.4	8.5	100353	Type 1 New
	Implante ZM4 MT	ZM44410M	RP	4.4	10	100353	Type 1 New
	Implante ZM4 MT	ZM44411M	RP	4.4	11.5	100353	Type 1 New
	Implante ZM4 MT	ZM44413M	RP	4.4	13	100353	Type 1 New
<b>ZIMMER DENTAL</b>	AdVent			4.5		100446	Type 35
	AdVent			5.7		100444	Type 34
	Eztetic Implant System CT31			3.1		100543	Type 60
	Screw-vent		3.5	3.3		100440	Type 32
	Screw-vent		3.5	3.7		100440	Type 32
	Screw-vent		4.5	4.7		100431	Type 27
	Spline			3.25		100436	Type 30
	Spline			3.75		100436	Type 30
	Spline			5.0		100436	Type 30
	SwissPlus			3.8		100425	Type 26
	SwissPlus			4.8		100350	Type 04
	Tapered Screw-vent		3.5	3.7		100440	Type 32
	Tapered Screw-vent		3.5	4.1		100440	Type 32
	Tapered Screw-vent		4.5	4.7		100431	Type 27
	Tapered Screw-vent		5.7	6.0		100442	Type 33
<b>ZUGA MEDICAL</b>	Bone Level			3.5		100554	Type 63
	Bone Level			5.0		100557	Type 65

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Angle	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>CAMLOG</b>	Comfour Bar abutment		Yellow		3.8		100579	Type 73
	Comfour Bar abutment		Purple		4.3		100579	Type 73
	Comfour Bar abutment		Blue		5.0		100580	Type 74
<b>COCHLEAR</b>	Baha BA 210					5.5	100487	Type 55
	Baha BA 210					8.5	100487	Type 55
	Baha BA 300, BIA 300					6.0	100487	Type 55
	Baha BA 300, BIA 300					9.0	100487	Type 55
	Baha BA 400, BIA 400					10.0	100487	Type 55
	Baha BA 400, BIA 400					12.0	100487	Type 55
	Baha BA 400, BIA 400					14.0	100487	Type 55
	Baha BA 400, BIA 400					6.0	100487	Type 55
	Baha BA 400, BIA 400					8.0	100487	Type 55
	Vistafix Standard abutment					3.0	100363	Type A3
	Vistafix Standard abutment					4.0	100363	Type A3
	Vistafix Standard abutment					5.5	100363	Type A3
	Vistafix Standard abutment					7.0	100363	Type A3
	Vistafix VXA 300					3.5	100363	Type A3
	Vistafix VXA 300					4.5	100363	Type A3
Vistafix VXA 300					6.0	100363	Type A3	
Vistafix VXA 300					7.5	100363	Type A3	
<b>NEOSS</b>	Access abutment						100483	Type 52
<b>NOBEL BIOCARE</b>	Brånemark system MUA		NP				100363	Type A3
	Brånemark system MUA		RP				100363	Type A3
	On1 Base CC		NP				100579	Type 73
	On1 Base CC		RP				100579	Type 73
	On1 Base CC		WP				100579	Type 73
	Mua		NP				100415	Type 25
	Mua		RP				100415	Type 25
Mua		WP				100415	Type 25	
<b>OTICON MEDICAL</b>	Ponto Implant abutment					4.0	100487	Type 55

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Angle	Diamètre (mm)	Longueur (mm)	N°Art.	SmartPeg
<b>STRAUMANN</b>	Coping for screw-retained abutment, Ti Bridge	023.2749	NC		3.5		100415	Type 25
	Coping for screw-retained abutment, Ti Bridge	023.4751	NC/RC		4.6		100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.2749	NC	17	4.6	2.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.2755	NC	17	4.6	4	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.0010	NC	17	4.6	5.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.2751	NC	30	4.6	2.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.2757	NC	30	4.6	4	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.0012	NC	30	4.6	5.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.4747	RC	17	4.6	2.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.4752	RC	17	4.6	4	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.4749	RC	30	4.6	2.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.4754	RC	30	4.6	4	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.0014	RC	17	4.6	5.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type A	022.0016	RC	30	4.6	5.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.2750	NC	17	4.6	2.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.2756	NC	17	4.6	4	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.0011	NC	17	4.6	5.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.2752	NC	30	4.6	2.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.2758	NC	30	4.6	4	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.0013	NC	30	4.6	5.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.4748	RC	17	4.6	2.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.4753	RC	17	4.6	4	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.4750	RC	30	4.6	2.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.4755	RC	30	4.6	4	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.0017	RC	30	4.6	5.5	100415	Type 25
	Screw-retained abutment, TAN, angled, Type B	022.0015	RC	17	4.6	5.5	100415	Type 25
	Screw-retained abutment, TAN, straight	022.2745	NC	0	3.5	1	100415	Type 25
	Screw-retained abutment, TAN, straight	022.2746	NC	0	3.5	2.5	100415	Type 25
	Screw-retained abutment, TAN, straight	022.2753	NC	0	3.5	4	100415	Type 25
	Screw-retained abutment, TAN, straight	022.2747	NC	0	4.6	1	100415	Type 25
	Screw-retained abutment, TAN, straight	022.2748	NC	0	4.6	2.5	100415	Type 25
	Screw-retained abutment, TAN, straight	022.2754	NC	0	4.6	4	100415	Type 25
	Screw-retained abutment, TAN, straight	022.4745	RC	0	4.6	1	100415	Type 25
	Screw-retained abutment, TAN, straight	022.4746	RC	0	4.6	2.5	100415	Type 25
	Screw-retained abutment, TAN, straight	022.4751	RC	0	4.6	4	100415	Type 25

Fournisseur d'implants	Modèle	N° Catalogue	Plateforme	Angle	Diamètre (mm)	Longueur (mm)	SmartPeg	N°Art.
<b>THOMMEN MEDICAL</b>	SPI System Bar abutment						100447	Type 36
	VARIOMulti abutment						100415	Type 25

\* Added the Implant Companies catalog no



### Un nouveau départ à chaque fois !

Afin d'éviter d'endommager l'implant, le matériau du SmartPeg est assez doux, ce qui limite la durée de vie du filetage. Ainsi, les SmartPegs sont à usage unique.

Rendez-vous sur [www.osstell.com/smartpegguide/](http://www.osstell.com/smartpegguide/) pour avoir le dernier Guide des références car celui-ci est régulièrement mis à jour.

